**Indications**

- Plantar fasciitis
- Heel spurs
- Achilles tendon injuries
- Ankle contracture

**Contraindications**

- Plantar flexion contracture greater than 10°
- Uncontrollable severe plantar flexion spasticity
- Extended weight bearing and ambulation

**Cleaning Instructions/Maintenance**

- When the shell, padding, or toe loop is slightly dirty, clean it using a damp cloth or a soft brush.
- Allow the shell, padding, or toe loop to air dry only.

**Instructions for use**

1. Remove both sides of the calf strap from the shell and remove the toe loop.
2. Slide foot into the night splint, seating the heel snugly in the heel opening.
3. Tighten both sides of the heel strap at the same time.
4. Tighten both sides of the toe strap at the same time.
5. Thread the calf strap. Tighten both sides of the calf strap at the same time.
6. Insert the big toe in the toe loop. The end of the toe should be aligned with the end of the toe loop.
7. Using the elastic ring, place the big toe in a neutral or slightly elevated position and attach the ring to one of the hooks. As treatment continues, the toe(s) may be elevated as comfort allows to permit a greater stretch to the plantar fascia.
8. Secure the toe loop strap to the bottom of the night splint.

**Caution**

Please always follow these instructions for use. In case of pain or swelling, discontinue use and consult your medical provider immediately. This device is designed to assist recovery and should be used exclusively under the instruction of and in consultation with a healthcare professional. DARCO International Inc. does not assume any responsibility for misuse of the product. Do not drive a vehicle while you are wearing this device. Use extreme caution when walking on all surfaces when using this device. Walk more slowly than normal. This device is for single-patient use only.

*To the patient/user: This device has been designed and manufactured by DARCO International, Inc. to be state of the art for treatment of a specific medical condition. Should any serious incident occur in relation to this device it should be reported to the manufacturer at the contact information listed on this document and the proper authority of the country or state where you are located. When treatment is complete, dispose of this device according to local laws and ordinances.*

**Manufactured for:**  
DARCO International, Inc.  
810 Memorial Blvd.  
Huntington, WV 25701, USA  
Phone 304.522.4883

**UK Responsible Person:**  
V-M Orthotics Ltd.  
Unit 25 Halesworth Business  
Centre Norwich Road  
Halesworth Suffolk IP19 8QJ  
United Kingdom




DARCO (Europe) GmbH  
Gewerbegebiet 18  
82399 Raisting, Germany  
Telefon +49 8807 9228 0  
www.darco.de

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