Important: This product is to be used **ONLY** under the supervision of a medical professional.

Use the PW series PegAssist™ with:

PW Series Off-loading Insole

PegAssist™

PegAssist™ Off-loading Insoles

The DARCO PegAssist™ System features a multi-indication removable-peg chassis that effectively off-loads the plantar aspect of the foot after surgery or when wounds are present.

- Removable pegs allow for localized off-loading of wounds and ulcerations of the foot.
- Pressure is reduced by as much as 60% while allowing patient to remain ambulatory.
- Each insole has a 15mm combination Plastazote® and multi-foam base with a 3mm PORON® Cushioning cover.
- Stabilizer board with low profile molded hook strip adheres to liner to prevent shifting while walking.

Contents:

- PegAssist™ Insole
- Stabilizer board
- Red marker

Video demonstration:


Proper assembly and application of the PegAssist™ Off-loading Insole

Scan QR code or visit:

darcointernational.com/pegassist

PW Series PegAssist™ Sizing:

Size:	Women:	Men:	Product Code:
XS	3 - 6	-	PW0
S	6.5 - 9	5 - 7.5	PW1
M	9.5 - 12	8 - 10.5	PW2
L	12.5+	11 - 13	PW3
XL	-	13.5+	PW4

Important: This product is to be used **ONLY** under the supervision of a medical professional.

Manufactured for:

DARCO International, Inc.
810 Memorial Blvd.
Huntington, WV 25701, USA
Toll-free: 800.999.8866
Phone: 304.522.4883
www.darcointernational.com

DARCO (Europe) GmbH
Gewerbegebiet 18
82399 Raisting, Germany
Tel. +49 (0) 88 07. 92 28-0
www.darco.de

Made in China

Instructions For Use

EN

Indications

- Off-loading, wound care management, pain management, foot trauma, post-op treatment

Contraindications

- Open fractures

Cleaning Instructions/Maintenance

- (Top layer only) When the PegAssist™ is slightly dirty, clean it with a damp cloth or a soft brush.
- Allow the PegAssist™ to air dry only.

Instructions for use

1. Separate walker straps and open the liner.
2. Insert the PW PegAssist™ Insole into the shoe with the **gray peg side facing down**. Ensure the insole is pushed all the way back and seated firmly against the footbed and heel.
3. Place a transparent barrier or dressing over the wound area. With the barrier in place, circle the wound or ulcer on the patient foot with the supplied marker.
4. Have the patient step onto the insole which will transfer the marking onto the top surface of the insole. (Remove barrier)
5. Remove the insole from the walker. Remove pegs beneath the marking to create the off-loading area. *Important: Remove the pegs slowly to avoid tearing the top layer of PORON®.*
6. Remove the adhesive backing from the stabilizer board that is included in the kit.
7. Align and attach the adhesive side of the stabilizer board to the **gray peg side** of the insole.
8. Insert the PW PegAssist™ Insole back into the walker at a downward angle. Ensure the stabilizer board is facing down, **blue side up**. Push all the way back so the heel of the insole meets the heel of the walker. Press the insole down to secure into the liner.
9. Insert foot and leg into walker and secure as per walker manufacturer instructions.

Caution

Please always follow these instructions for use. In case of pain or swelling, discontinue use and consult your medical provider immediately. This device is designed to assist recovery and should be used exclusively under the instruction of and in consultation with a healthcare professional. DARCO International Inc. does not assume any responsibility for misuse of the product. It is advisable to wear dressings or socks with the PegAssist™. Do not drive a vehicle while you are wearing this device. Use extreme caution when walking on all surfaces when using this device. Walk more slowly than normal. This device is for single-patient use only.

To the patient/user: This device has been designed and manufactured by DARCO International, Inc. to be state of the art for treatment of a specific medical condition. Should any serious incident occur in relation to this device it should be reported to the manufacturer at the contact information listed on this document and the proper authority of the country or state where you are located. When treatment is complete, dispose of this device according to local laws and ordinances.

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