Body Armor® PFA Walker

MD	EN Medical Device DE Medizinprodukt ES Producto sanitario FR Dispositif médical IT Dispositivo medico NL Medisch hulpmiddel DA Medicinsk udstyr SV Medicinteknisk produkt PL Urządzenie medyczne SL Medicinski pripomoček HR Medicinski proizvod PT Dispositivo médico	UDI	EN Unique device identifier DE Eindeutige Produktidentifizierung ES Identificación única del producto FR Identifiant unique des dispositifs IT Identificazione univoca del dispositivo NL Unieke identificatie van het hulpmiddel DA Unik produktidentificering SV Entydigt produkt-ID PL Unikalny identyfikator produktu SL Edinstvena identifikacija pripomočka HR Jedinstveni identifikacijski proizvod PT Identificador único do dispositivo
-----------	---	------------	---

EN	Manufacturer	EN	Single patient multiple use
DE	Hersteller	DE	Mehrfachanwendung an einem einzelnen Patienten
ES	Fabricante	ES	Reutilizable en un mismo paciente
FR	Fabricant	FR	Utilisation multiple sur un seul patient
IT	Fabbricante	IT	Applicazione multipla su un solo paziente
NL	Fabrikant	NL	Meerdere keren te gebruiken bij één patient
DA	Producent	DA	Fleere anvendelser på en enkelt patient
SV	Tillverkare	SV	Användning flera gånger på samma patient
PL	Producent	PL	Wielokrotne użycie u jednego pacjenta
SL	Proizvajalec	SL	Večkratna uporaba pri enem bolniku
HR	Proizvođač	HR	Višestruka primjena na jednom pacijentu
PT	Fabricante	PT	Aplicações múltiplas num único doente

EN	Instructions for use
Body Armor® PFA Walker Foot Stump Orthosis	
Important notes If there is a serious deterioration in your state of health as a result of using the medical appliance, please report this to your specialist dealer or the manufacturer as soon as possible. If you report this to the specialist dealer, the dealer will forward this immediately. The reports are agreed with and processed by the competent authorities.	
Please always follow these instructions for use. In case of adverse effects (e.g. skin swelling, pressure points, disorders of blood circulation, etc.) please tell your doctor. Function and optimum protection are no longer guaranteed if the product is handled incorrectly. This device assists recovery and should be used exclusively under the instruction of and in consultation with your doctor or a health-care professional. DARCO (Europe) GmbH does not assume any responsibility. Do not wear the product in direct contact with injured, irritated, damaged or broken skin. This device is not recommended for direct skin contact. You should use bandages, dressings or socks.	
You must not drive a vehicle while you are wearing this product. Use extreme caution when climbing stairs (e.g. also stepping up onto curbs) and walking on uneven surfaces. Be particularly cautious when walking in wet, slippery and snowy conditions. Walk more slowly than usual.	
Note: This device is allowed for single-patient use only. It does not need to be cleaned or reprocessed between uses.	
Purpose	Foot Stump Orthosis as an interim treatment for stabilisation and offloading with simultaneous mobilisation of the patient.
Indications / Contraindications	<p>Indications Postoperative for partial foot amputations (e.g. Chopart or Lisfranc amputations) due to</p> <ul style="list-style-type: none"> > Diabetic foot syndrome > PAOD (peripheral arterial occlusive disease) > Trauma <p>Contraindications No contraindications are known.</p>
Design features	Stable, ankle-spanning shaft with neutral insole for offloading of the stump. The integrated and adjustable air cushions enable individual adjustment. The product contains a flexible ventral support strap. The variably adjustable hook and loop fasteners with a guiding tab provide good circular stability.
Combination with other products	No specifications
Maintenance	The product does not require maintenance.

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EN Use

- Photo 1:** Open the hook and loop fasteners on the orthosis and remove the grey inner shoe. Open the hook and loop fasteners of the inner shoe and place this around the leg. Fasten the inner shoe using the hook and loop fasteners. First close the hook and loop tape over the instep, then the lower hook and loop strap as numbered in the photo.
- Photo 2:** Now place the insole into the walker. To insert the foot securely into the walker, spread the medial and lateral parts of the shell with both hands. Step onto the insole with the grey inner shoe. The insole can be adjusted individually by qualified staff.
- Photo 3:** After inserting the foot, pass the straps into the guiding tabs and adjust them to the desired stability.
- Photo 4:** Now pump air through the round pump on the front of the walker into the air chambers of the grey inner shoe (A). The air is distributed evenly. The compression achieves additional stability of the ankle joint. To reduce the volume of the air chambers, release the air using the silver valve (B).

Note: If impaired sensation, swelling or pain occurs, the pressure should be reduced by releasing air from the chambers. Should persistent pain occur, the treatment should be interrupted and a doctor should be consulted.

Technical data

3 sizes
M = stump length up to 19.5 cm
L = stump length up to 21 cm
XL = stump length up to 22.5 cm
Universal fit for right and left.

Colour
Grey-Blue

Material
Outsole walking surface: EVA
Plastic shell: Polypropylene, PU
Fiber liner: Polyester
Buckles: Iron
Rivets: Iron
Straps: Nylon / PVC
Insole: EVA and MTPU

Working life / Lifespan of the product

The lifespan of the medical device is determined by natural wear and tear due to proper handling/cleaning. The medical device is designed to be used for the usual duration of use prescribed by a doctor without showing a reduction in safety or performance. For safety, the functionality of the product should be checked every 6-8 weeks together with a medical professional or at a specialist dealer.

Care instructions

If the product is slightly soiled, you can wipe it with a damp cloth or clean it with a soft brush. To clean the shell, use lukewarm water (30° is possible) and mild soap or mild detergent if necessary. Rinse and dry the liner promptly. Do not dry this in a dryer or on a radiator. Let it dry in the air at room temperature.

Disposal

Dispose of after use. Please dispose of the packaging in an environmentally friendly manner. Dispose of this product after use in accordance with applicable local laws and regulations.

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