

The Body Armor® Sport New By DARCO

Perfect stability along with controlled mobility!



The Body Armor® Sport is an ankle joint orthosis which is applied during the acute phase of injury as well as during the process of rehabilitation.



DARCO



1. Elements of stability

U-shaped heel plate is anatomically shaped to the medial /lateral sides of the foot.

2. Orthosis sole

Additional space in the shoe due to shortened orthosis sole.

3. DARCO Y-closure system

Straps may be cut to desired length and resealed with the Y-closure system.

4. Talus shift strap

Protection of the anterior talo-fibular ligament by means of the the talus shift strap from plantar /lateral to dorsal /medial. (Prevents the frontal shift of the talus.)

5. Parallel hinge joints

Physiologically correct movement of the foot by parallel and thus anatomically correct positioning of the hinge joints.

6. Customizable liner

Can be modified to transition from acute to supportive brace.

Features and Benefits:

- > Post-traumatic and preventive care.
- > Anatomically shaped stirrups ensures a perfect fit.
- > Customizable.
- > Physiological correct movement of the foot.
- > Prevents the frontal shift of the talus.

Indications

- > Distorsion trauma of the upper ankle joint
- > Insufficiency of the capsules and ligaments of the upper ankle joint
- > Prevention of chronic instability

Sizing Chart

Size	Item-No.	Men's	Womens
Small	BAS1L	Below 7	5-8
	BAS1R		
Medium	BAS2L	7.5-10.5	9-12
	BAS2R		
Large	BAS3L	11-14	13+
	BAS3R		



DARCO international / 810 Memorial Blvd. / Huntington, WV 25701
 Toll Free 800-999-8866 / Toll Free Fax 877-201-7082 / Voice 304-522-4883
 info@darcointernational.com / www.darcointernational.com