

# FX Pro™ Walker

Ultra Low Profile, Lightweight Walker



The re-designed FX Pro™ Walker features a more rigid shell with a soft direct injection outsole for better traction and less pressure at heel strike. Our low profile design eliminated leg length discrepancies associated with other walkers and is built with a wider platform to better accommodate edema or large bandages. The new “soft touch” band edging the stirrups and softer strapping materials provides greater patient comfort and compliance.





## FX Pro™ Features and Benefits:

- › **Low Profile**, eliminates height difference between walker and normal shoe without altering gait and reducing hip and knee pain.
- › **Padded liner** is ventilated which promotes air circulation and aids in keeping the lower leg cool and comfortable.
- › **Anatomically contoured “spatula” stirrups** immobilize the ankle comfortably.
- › **Pressure moldable, shock absorbing insole** conforms to the foot for maximum comfort.

## FX Pro™ Indications:

- › Post operative immobilization
- › Strains and sprains of the lower leg, ankle or foot
- › Off-loading wounds or ulcerations when used with PegAssist™ Insole (PW Series)

## Sizes and Product Codes



High Version:	Women:	Men:	Product Code
XS	3.5 - 6	-	FX0
S	6.5 - 9	5 - 7.5	FX1
M	9.5 - 12	8 - 10.5	FX2
L	12.5+	11 - 13	FX3
XL	-	13.5+	FX4



Low Version:	Women:	Men:	Product Code
XS	3.5 - 6	-	FXS0
S	6.5 - 9	5 - 7.5	FXS1
M	9.5 - 12	8 - 10.5	FXS2
L	12.5+	11 - 13	FXS3
XL	-	13.5+	FXS4

Available through:



Headquarters:

DARCO International | 810 Memorial Blvd. | Huntington, WV 25701 | USA

Toll Free voice 800-999-8866 | Fax 877-201-7082 | Voice 304-522-4883 | Fax 304-522-0037

info@darcointernational.com | www.darcointernational.com

*As with all DARCO products, a portion of the profits from the sale of these items will be donated to encourage and support continuing education in our profession.*