


# MedSurg™ Shoe

## Post Operative Protective Footwear



The DARCO MedSurg™ Shoe is the most versatile post operative shoe on the market today. It's fully customizable insole is double-padded and removable to accommodate the PegAssist™ off-loading insole.

**DARCO**



## MedSurg™ Features and Benefits:

- › Square Toe Design acts as a bumper to provide extra protection and facilitates better left/right fit.
- › MetaShank™ protection provides more rigid control under the metatarsal heads.
- › Stapless Closure eliminates buckle pressure on the forefoot and can expand to accommodate dressings up to 17 inches in circumference.
- › Rocker Sole reduces plantar pressure on the heel and forefoot by over 25%.
- › Fully compatible with DARCO's PegAssist™ off-loading insole.

## MedSurg™ Indications:

- › After foot surgery or when K-wires are present
- › When large bandages or dressings are required
- › As a platform for the PegAssist™ insole when off-loading wounds

Women's Item Number	Men's Item Number	Description	Size	Color	Women's Size	Men's Size
MQW1B	MQM1B	MedSurg	Small	Black	4 - 6	6 - 8
MQW2B	MQM2B	MedSurg	Medium	Black	6.5 - 8	8.5 - 10
MQW3B	MQM3B	MedSurg	Large	Black	8.5 - 10	10.5 - 12
	MQM4B	MedSurg	X-Large	Black		12.5 - 14

Available through:



Headquarters:

DARCO International | 810 Memorial Blvd. | Huntington, WV 25701 | USA

Toll Free voice 800-999-8866 | Fax 877-201-7082 | Voice 304-522-4883 | Fax 304-522-0037

info@darcointernational.com | www.darcointernational.com

*As with all DARCO products, a portion of the profits from the sale of these items will be donated to encourage and support continuing education in our profession.*