Body Armor® Cast Shoe

Cast Boot



The redesigned DARCO Body Armor® Cast Shoe features improvements in function and design. It is specifically engineered to shield the cast from moisture, shock and everyday wear. The shoe can be used with all casting materials and bulky bandages while sporting a fresh, modern design that will appeal to all patients.









Body Armor® Cast Shoe Features and Benefits:


- > The improved design features sturdier low profile buttons, a more secure locking mechanism, a stronger bungee system, and a redesigned body that better adapts to the foot.
- > EVA Unibody Design provides unparalleled shock absorption while protecting the cast from the elements.
- > Unique Bungee Closure System is easy to use, easy to adjust and effectively secures the cast in the shoe.
- > Built in Toe Guard protects the toes and prevents debris from entering the toe area
- > Style that sets it apart; Function that puts it in a class of its own.
- > Six Sizes to accommodate any foot, while reducing inventory.
- > Suggested Code: L3260

Body Armor® Cast Shoe Indications:

- Cast protection
- Bulky bandages protection

Body Armor® Cast Shoe Sizing in Inches:

Size:	Length (A):	Width (B):	Width (C):	Product Codes
XS	9.4"	3.9"	3.5"	BACS0
S	10"	4.4"	4.1"	BACS1
M	10.6"	4.7"	4.3"	BACS2
L	11.4"	5.1"	4.6"	BACS3
XL	12.6"	5.5"	4.8"	BACS4
2XL*	13.7"	5.9"	5.1"	BACS5



*2XL coming soon.

Available through:

