Body Armor[®] Night Splint Dorsal Night Splint



The DARCO Body Armor[®] Night Splint is the only night splint that actively engages the windlass mechanism of the foot, resulting in a specific and sustained stretch to the plantar fascia and a sustained low load stretch to the flexor tendons, Achilles tendon and calf muscles.





Body Armor® Night Splint Features and Benefits:

- Toe Strap dorsiflexes the 1st MPJ to provide a specific and sustained stretch to the Plantar Fascia and a sustained low load stretch to the flexor tendons, Achilles tendon and calf muscles.
- > Loading the Windlass Mechanism naturally brings the foot/ankle to a near 90 degree angle thus eliminating the "lever effect" and the resulting pressure on the apex of the ankle present in other dorsal models.
- Elastic Connector Ring allows for natural flexion of the foot and toes during sleep and gently returns the foot to the dorsiflexed position as muscles relax.
- > Low Profile Design allows for maximum user comfort and compliance.
- > One-size-fits-all means reduced inventory and ensures that the right size is always on hand.

Body Armor® Night Splint Indications:

 Plantar Fasciitis 	 Achilles Tendonitis 	 Metatarsalgia
› Heel Spurs	 Achilles Tendinosis 	 Ankle Contracture

Item Number	Description	Size	Color
BADS	Dorsal Night Splint	Universal	White

Available through:



Headquarters:

DARCO International | 810 Memorial Blvd. | Huntington, WV 25701 | USA Toll Free voice 800-999-8866 | Fax 877-201-7082 | Voice 304-522-4883 | Fax 304-522-0037 info@darcointernational.com | www.darcointernational.com

As with all DARCO products, a portion of the profits from the sale of these items will be donated to encourage and support continuing education in our profession.