
Body Armor® Toe Guard

A cost-effective solution.



The DARCO Body Armor[®] Toe Guard provides extra protection for the toe region during the recovery process. It works particularly well when k-wires are present or anytime when accidental contact could cause pain or be detrimental to the healing process.





Body Armor® Toe Guard Features and Benefits:

- > Two models are available. The Round Toe Guard fits all round toe post operative shoes. The Square Toe Guard fits Darco square toe shoes and walkers.
- > Hard ABS Acrylonitrile Butadiene Styrene plastic withstands blows and strikes to the toe region in everyday situations.
- > Moldable side wings can be bent or molded to conform to any foot size.
- > Adhesive backing stays strong and secure in shoes for patient comfort and security.

Body Armor® Toe Guard Indications:

- Post Operative Healing
- Presence of K-Wires

Trauma Care

Diabetic Foot/Neuropathy

Toe Guard Sizing Chart:

Square	Round	Size
BATGQ1	BATGS1	Women's Small Women's Medium Women's Large
BATGQ2	BATGS2	Men's Small Men's Medium Men's Large Men's X-Large

Available through:



Headquarters:

DARCO International | 810 Memorial Blvd. | Huntington, WV 25701 | USA Toll Free voice 800-999-8866 | Fax 877-201-7082 | Voice 304-522-4883 | Fax 304-522-0037 info@darcointernational.com | www.darcointernational.com

As with all DARCO products, a portion of the profits from the sale of these items will be donated to encourage and support continuing education in our profession.