

WCS™ Closed-Toe

Wound Care Shoe System



The WCS™ Closed-Toe Wound Care Shoe System is the product of choice when dealing with Wagner Grade 3 or greater ulcerations. The four multi-density insoles can be modified, effectively decreasing pressure from the plantar aspect of the foot.



DARCO



Replacement insoles available. Sold separately.

WCS™ Closed-Toe Features and Benefits:

- › **Four Multi-density Insoles** allow for the ultimate in customization and precisely targeted off-loading.
- › **Ultra-soft Plastazote® Lining** provides a virtually friction-free contact surface that won't exacerbate existing wounds.
- › **Sold in Pairs** to eliminate discrepancies in shoe height and the potential associated knee and hip pain.

Insoles Cross Section

- Plastazote®
Thickness 1/4" Durometer 26 →
- Medium Density EVA
Thickness 1/4" Durometer 42 →
- Firm Density EVA
Thickness 1/4" Durometer 50 →
- PORON® Cushioning
Thickness 1/8" Durometer 29 →



WCS™ Closed-Toe:

- › Wagner Grade 3 or greater ulcerations
- › Bony prominences
- › Long-term treatment of wounds and ulcerations
- › Maximum off-loading while still allowing for patient ambulation

Item Number	Size	Color	Women's Size	Men's Size
WCC0B	X-Small	Black	4 - 5.5	-
WCC1B	Small	Black	6 - 8.5	6 - 7.5
WCC7B	Small/Medium	Black	9 - 9.5	8 - 8.5
WCC2B	Medium	Black	10 - 10.5	9 - 9.5
WCC3B	Large	Black	11 - 12.5	10 - 11.5
WCC4B	X-Large	Black	13 - 14.5	12 - 13.5
WCC5B	2X-Large	Black	15 - 17	14 - 16

Available through:



Headquarters:

DARCO International | 810 Memorial Blvd. | Huntington, WV 25701 | USA

Toll Free voice 800-999-8866 | Fax 877-201-7082 | Voice 304-522-4883 | Fax 304-522-0037

info@darcointernational.com | www.darcointernational.com