


HeelWedge™ Healing Shoe


Effective Off-Loading Footwear



Standard
Surgical Shoe



DARCO
HeelWedge™



The DARCO HeelWedge™ is clinically proven to off-load pressure from the heel by shifting weight to the mid and forefoot to promote faster healing after surgery, trauma or when wounds or ulcerations are present on the heel.

DARCO




HeelWedge™ Features and Benefits:

- › Reduces Heel pressure by **over 26%**!
- › **Square Toe Design** acts as a bumper to protect the toes and provides better universal Left/Right fit.
- › **Ankle Strap** seats foot firmly in shoe to reduce heel slippage and eliminate friction.
- › **Removable Forefoot Closure** provides security while eliminating buckle pressure.
- › **Zoned Outsole** has an aggressive tread under the mid-foot where traction is needed most.
- › **Removable Insole** features twice the padding of standard insoles and can be modified as needed. Insole can be totally removed to accommodate DARCO's innovative **PegAssist™** customizable off-loading insole for an even more targeted off-load.

HeelWedge™ Indications:

- › Rearfoot trauma
- › Wounds or ulcerations present on the heel area.
- › Post Surgical healing for either soft tissue or bony structure of the heel

Item Number	Description	Size	Color	Women's Size	Men's Size	Optional PegAssist™ Insole
HQ0B	HeelWedge	X-Small	Black	3 - 5.5	-	PQ0
HQ1B	HeelWedge	Small	Black	6 - 7.5	5 - 6.5	PQ1
HQ2B	HeelWedge	Medium	Black	8 - 9	7 - 8	PQ2
HQ3B	HeelWedge	Large	Black	9.5 - 11	8.5 - 10	PQ3
HQ4B	HeelWedge	X-Large	Black	11.5 - 13	10.5 - 12.5	PQ4



Pair with the PQ PegAssist™ to enhance offloading.

Available through:



Headquarters:
 DARCO International | 810 Memorial Blvd. | Huntington, WV 25701 | USA
 Toll Free voice 800-999-8866 | Voice 304-522-4883
 darco@darcointernational.com | www.darcointernational.com