



Body Armor[®]

Stirrup Walker

Engineered for Comfort

Featuring DARCO's patented PowerPods[®] Outsole, the Body Armor[®] Stirrup Walker delivers on comfort throughout the healing process.

Made from a continuous layer of high-rebounding EVA foam, DARCO PowerPods[®] hold their shape and thickness over time to deliver lasting comfort and protection you can rely on.

www.darcointernational.com



Featuring PowerPods® Outsole

DARCO's PowerPods® sole features a continuous layer of EVA foam, a material commonly used in athletic footwear. This proprietary* sole construction improves the performance and durability of the footwear to provide comfort that lasts.



Bilateral



For diabetics



Time saving

Indications

- Postoperative use
- Trauma care
- Stress fractures
- Soft tissue injury

Contraindications

- Open fractures
- Poor circulation


Product Options

- Air, Non-Air
- Short, Tall
- 5 sizes

Suggested Code:
L4386 (non-pneumatic)
L4360 (pneumatic)



Sizes + Product Codes



Size	Women	Men				Tall		Short	
			A	B	C	Air	Non-Air	Air	Non-Air
XS	3.5 – 6.0	–	8.7	4.0	3.1	BPT0G	BST0G	BPS0G	BSS0G
S	6.5 – 9.0	5.0 – 7.5	9.1	4.1	3.3	BPT1G	BST1G	BPS1G	BSS1G
M	9.5 – 12.0	8.0 – 10.5	11.0	4.4	3.5	BPT2G	BST2G	BPS2G	BSS2G
L	12.5+	11.0 – 13.0	11.2	4.6	3.7	BPT3G	BST3G	BPS3G	BSS3G
XL	–	13.5+	11.7	4.8	4.1	BPT4G	BST4G	BPS4G	BSS4G