



New
Product!



DUO™ Pro

Wound Care Shoe

Experience superior wound care with our lightweight & durable **all-in-one solution.**

Designed with state-of-the-art materials and an innovative off-loading insole, DUO™ Pro is 27% lighter and reduces pressure by 40% compared to competitor post-op shoes. Built to last 4x longer, this comprehensive solution includes a specialized shoe, customizable insole, and wound marker for custom patient care.

www.darcointernational.com





Featuring Aircore INSOLE Technology

DARCO's Aircore™ Insole further enhances personalized patient care by enabling the creation of a custom off-loading zone.



Bilateral



For diabetics



Time saving

Indications

- Treatment of diabetic foot syndrome
- To facilitate healing of open wounds and minor ulcerations
- Pressure off-loading
- Up to Wagner grade 3 ulcerations
- Bony prominences

Contraindications

- Open fractures
- Ulcerations greater than Wagner grade 3

Suggested Code: L3260

The Benefits

Dual Hardness Technology

Features a traction outsole layer for enhanced grip and stability, paired with a cushioning outsole layer designed to provide shock absorption and reduce pressure.

Pressure Reduction

Two layers of specially engineered foaming EVA offer superior cushioning, shock-absorption, support and stability—with up to 40% more pressure reduction than competitor post-op shoes. The Aircore™ Insole further enhances personalized patient care by enabling the creation of a custom off-loading zone.

Longevity

Designed to last four times longer than competitor post-op shoes, providing unmatched durability and value.



Sizes + Product Codes

Size	Women	Men	Insole			SKU #
			A	B	C	
XS	3.5 – 6.0	–	9.4	3.9	3.0	DP0B
S	6.5 – 7.5	5.0 – 6.5	10.2	4.3	3.3	DP1B
M	8 - 9	7 - 8	10.8	4.4	3.4	DP2B
L	9.5 - 11	8.5 - 10	11.8	4.6	3.7	DP3B
XL	11.5 - 13	10.5 - 12.5	12.3	4.9	3.9	DP4B

