



The DARCO **MedSurg™** now features a durable EVA outsole that is 20% lighter than previous models, with a specially designed zoned outsole providing improved traction and 150% better rebound to retain its shape and thickness over time. It is the most versatile post-op shoe on the market today.

**DARCO**



**MedSurg™**  
Post-Op Shoe



Bilateral



For diabetics



Time saving


## Indications

- Postoperative for stabilization and off-loading
- Postoperative healing (e.g. foot osteotomies)
- Trauma care
- Presence of K-wire

## Contraindications

- Open fractures

**Suggested Code: L3260**



## The Benefits

### Zoned Outsole

The DARCO MedSurg™ features an EVA outsole 20% lighter than previous models, with our specially designed zoned outsole providing improved traction.

### Highly Durable

Durable EVA outsole withstands compression with 150%\* better rebound to retain shape and thickness over time.

### Offloading Capabilities

The MedSurg™ post op shoe with bamboo Metashank™ provides added support under the forefoot while walking, and reinforces the outsole increasing durability. Combine with the PTQ PegAssist™ offloading insole to transform our post-op shoe into a diabetic offloading solution.



## Sizes + Product Codes

Size	Women	Men	Insole (Women)			Insole (Men)			Women	Men
			A	B	C	A	B	C	SKU #	SKU #
<b>S</b>	4 - 6	6 - 8	9.2	3.5	2.2	10.3	4.1	2.5	MQW1B	MQM1B
<b>M</b>	6.5 - 8	8.5 - 10	9.6	3.7	2.2	10.8	4.3	2.6	MQW2B	MQM2B
<b>L</b>	8.5 - 10	10.5 - 12	8.2	3.8	2.3	11.4	4.5	2.8	MQW3B	MQM3B
<b>XL</b>	-	12.5 - 14	-	-	-	12	4.9	3.0	MQW4B	MQM4B
<b>Pediatric</b>		<b>Shoe Size: 12 -1</b>	-	-	-	-	-	-	MQPB	-