

DARCO



OrthoWedge™

Healing Shoe

The DARCO OrthoWedge is clinically proven to reduce weight bearing pressure on the forefoot which promotes faster healing after surgery, trauma or when forefoot wounds or ulcerations are present.

www.darcointernational.com





Compatible with PQ PegAssist Insoles
 Pair with the PQ PegAssist to enhance offloading

The Benefits

- Reduces forefoot pressure by as much as 57%!
- Square-toe design acts as a bumper to protect the toes.
- Ankle strap seats foot firmly in shoe to reduce heel slippage and eliminate friction.
- Removable forefoot closure eliminates buckle pressure.
- Zoned outsole has an aggressive tread under the midfoot.
- EVA insole features twice the padding of standard insoles and can be modified as needed.
- Optional PegAssist Insole can be added to provide enhanced offloading capabilities when treating diabetic wounds.

Indications

- Post surgical healing
- Forefoot trauma
- Wounds or ulcerations present under metatarsal heads and toes
- Any situation that requires a pressure off-load of the forefoot

Contraindications

- Open fractures

This product is not designed to be worn in pairs.

Suggested Code: L3260



Sizes & Product Codes

OQ™ Size	Women's Size	Men's Size	Item Number	Optional PegAssist™ Insole
X-Small	3 - 5.5	-	OQ0B	PQ0
Small	6 - 7.5	5 - 6.5	OQ1B	PQ1
Medium	8 - 9	7 - 8	OQ2B	PQ2
Large	9.5 - 11	8.5 - 10	OQ3B	PQ3
X-Large	11.5 - 13	10.5 - 12.5	OQ4B	PQ4

Version 1_1-2026