

Body Armor[®] PFA Walker VS ProTerm



Highlights



Larger Footbed Profile

Footbed profile is longer and wider, making it better-suited for recovery from TMA procedures



More Comfortable

PowerPods construction reduces amount of materials used, resulting in a lighter walker that reduces leg length discrepancy



Accommodates Heavy Bandaging

Walker shaft circumference is larger to better accommodate heavy bandaging and edema

Product Information

Pricing:

- DARCO Direct Price : \$175

Suggested Code:

- Pneumatic: L4360*, ~\$250

Available Through:

- Distributors
- DARCO Direct

	PFA Walker	ProTerm	
	Size: Medium	Size: Medium	
Floor to Top of Insole	50 mm	60 mm	16% lower profile
Insole Length	200 mm	165 mm	21% longer footbed
Insole Width at Forefoot	132 mm	95 mm	39% wider forefoot opening
Shaft Circumference	470 mm	410 mm	Accommodates 14% larger leg circumference
Total Weight	687 g	845 g	18% lighter

*Application Pending