



Wound Care Shoe System™

Open-Toe and Closed Toe

The WCS™ Open-Toe Wound Care Shoe System is the product of choice when dealing with lesions through Wagner Grade 3 ulcerations. The four multi-density insoles can be modified, effectively decreasing pressure from the plantar aspect of the foot.


www.darcointernational.com





Wound Care Shoe System™ Insoles

Replacement insoles available. Sold separately.



Indications

- Up to Wagner Grade 3 ulcerations
- Bony prominences
- Long-term treatment of wounds and ulcerations
- Maximum off-loading while still allowing for patient ambulation

Contraindications


- Open fractures
- Ulcerations greater than Wagner grade 3

Suggested Code: A9270

The Benefits


- Four Multi-density Insoles allow for the ultimate in customization and precisely targeted off-loading.
- Ultra-soft Plastazote® Lining provides a virtually friction-free contact surface that won't exacerbate existing wounds.
- Sold in Pairs to eliminate discrepancies in shoe height and the potential associated knee and hip pain.

Insoles Cross Section



Sizes & Product Codes

WCS™ Size	Women's Size	Men's Size	Item Number Open Toe	Item Number Closed Toe
X-Small	4 - 5.5	-	WCS0B	WCC0B
Small	6 - 8.5	6 - 7.5	WCS1B	WCC1B
Small/Medium	9 - 9.5	8 - 8.5	WCS7B	WCC7B
Medium	10 - 10.5	9 - 9.5	WCS2B	WCC2B
Large	11 - 12.5	10 - 11.5	WCS3B	WCC3B
X-Large	13 - 14.5	12 - 13.5	WCS4B	WCC4B
2X-Large	15 - 17	14 - 16	WCS5B	WCC5B



Version 1_1-2026

