

DARCO



HeelWedge™

Healing/Off-loading Shoe

The DARCO HeelWedge™ is clinically proven to off-load pressure from the heel by shifting weight to the mid and forefoot to promote faster healing after surgery, trauma or when wounds or ulcerations are present on the heel.

www.darcointernational.com





Compatible with PQ PegAssist Insoles
 Pair with the PQ PegAssist to enhance offloading

The Benefits

- Reduces heel pressure by over 26%!
- Square-toe design acts as a bumper to protect the toes.
- Ankle strap seats foot firmly in shoe.
- Removable forefoot closure eliminates buckle pressure.
- Zoned outsole has an aggressive tread under the midfoot where traction is needed most.
- EVA insole features twice the padding of standard insoles and can be modified.
- Optional PegAssist Insole can be added to provide enhanced offloading capabilities when treating diabetic wounds.

Indications

- Post surgical healing
- Wounds or ulcerations
- Any situation that requires a pressure off-load of the rearfoot
- Rearfoot trauma

Contraindications

- Open fractures

This product is not designed to be worn in pairs.

Suggested Code: L3260



Sizes & Product Codes

| OQ™ Size | Women's Size | Men's Size | Item Number | Optional PegAssist™ Insole |
|----------|--------------|-------------|-------------|----------------------------|
| X-Small | 3 - 5.5 | - | HQ0B | PQ0 |
| Small | 6 - 7.5 | 5 - 6.5 | HQ1B | PQ1 |
| Medium | 8 - 9 | 7 - 8 | HW2B | PQ2 |
| Large | 9.5 - 11 | 8.5 - 10 | HQ3B | PQ3 |
| X-Large | 11.5 - 13 | 10.5 - 12.5 | HQ4B | PQ4 |

Version 1_3_2026