

New



AllRound Shoe™

Closed Healing Shoe

The AllRound Shoe™ is an ideal choice for patients seeking a discreet, closed-toe healing solution that can be worn as a single or pair. Designed with a rounded toe box, it eliminates unnecessary bulk while still delivering maximum protection. The closed-toe construction shields the foot from moisture and debris, making it a great option for patients ready to return to work while continuing to heal responsibly. A fully accessible wide opening allows for easy entry—even over bulky dressings—while the seamless interior helps prevent skin irritation and breakdown. Finished with an EVA outsole for reliable shock absorption and slip resistance, the AllRound Shoe™ offers less bulk, more protection, and everyday comfort.

DARCO

The Benefits

Seamless Interior

- Constructed without seams or ridges to minimize friction and prevent harm to sensitive skin.

Easy Entry

- Fully accessible entry with wide opening to benefit patients who struggle with mobility
- Up to 45cm (17.7 inches) of midfoot accommodation for patients with swelling, heavy bandaging or custom orthotics—36% larger than standard post-operative solutions

Less Bulk. More Protection.

- Sleek, round toe box eliminates extra bulk & reduces trip hazards
- MetaShank for added stabilization of the plantar surface of the foot
- EVA-outsole for shock absorption + anti-slip
- Closed toe design to protect against moisture and debris



Bilateral



No Shoe Balancer
Required



Lightweight



Wide opening for easy foot insertion



Anti-slip outsole

Indications

- For postoperative care
- For wound care
- Arthritis

Contraindications

- No contraindications are known

Suggested Code: L3260

Adult Sizes + Product Codes

	US Women's	US Men's	in inches			SKU #
			A	B	C	
XS	Up to 7.5	Up to 6	9.8	3.4	2.8	ARS0
S	8 - 9	6.5 - 7.5	10.3	3.6	2.9	ARS1
M	9.5 - 10.5	8 - 9	10.7	3.7	3.0	ARS2
L	11 - 12	9.5 - 10.5	11.1	3.9	3.2	ARS3
XL	12.5 - 13.5	11 - 12	11.7	4.1	3.3	ARS4
2XL	14 - 14.5	12.5 - 13	12.4	4.3	3.5	ARS5

