



Body Armor®

Toe Guard

The DARCO Body Armor® Toe Guard provides extra protection for the toe region during the recovery process. It works particularly well when k-wires are present or anytime when accidental contact could cause pain or be detrimental to the healing process.



The Benefits

- **Two models available.** The Round Toe Guard fits all DARCO round toe post operative shoes. The Square Toe Guard fits DARCO square toe shoes and walkers.
- **Hard ABS** - Acrylonitrile Butadiene Styrene plastic withstands blows and strikes to the toe region in everyday situations.
- **Moldable side wings** can be bent or molded to conform to any foot size.
- **Adhesive backing** stays strong and secure in shoes for patient comfort and security



Sizes + Product Codes

Shoe Size	Square	Round
Women's Small - Large	BATGQ1	BATGS1
Men's Small - X LargeS	BATGQ2	BATGS2

Item sold as one (1) guard. Suggested Code: L9900

Indications

- Post Operative Healing
- Trauma Care
- Presence of K-Wires
- Diabetic Foot/Neuropathy