



A night splint that works with the foot—not against it. The Hallux Sling™ dorsiflexes the great toe to stretch the plantar fascia, naturally positioning the foot at 90° while reducing ankle pressure. Its lightweight, breathable design improves comfort, leading to better adherence and more effective relief for patients.

DARCO



Body Armor®

Dorsal Night Splint

The Benefits

- **Lightweight, breathable design** for comfortable use while at rest
- **Engaging the Windlass Mechanism** naturally brings the foot/ankle to a near 90 degree angle thus eliminating the “lever effect” and the resulting pressure on the apex of the ankle present in other dorsal models.
- **Elastic Connector Ring** allows for natural flexion of the foot and toes during sleep and gently returns the foot to the dorsiflexed position as muscles relax.
- **Low Profile Design** allows for maximum user comfort and compliance.
- **One-size-fits-most** means reduced inventory and ensures that the right size is always on hand.

Item sold as one (1) night splint. Fits right or left foot.
Universal sizing.

Sizes + Product Codes

Item #	Size	Color
BADS	Universal	White



Elastic connector ring allows you to find the right stretch for you

Indications

- Plantar Fasciitis
- Heel Spurs
- Achilles Tendon Injuries
- Ankle Contracture



Bilateral

Suggested Code: L4396



Naturally brings the foot to a 90-degree angle