



Body Armor®

Sport Ankle Brace

The Body Armor® Sport Ankle Brace is an ankle joint orthosis which is applied during the acute phase of injury as well as during the process of rehabilitation.



The Benefits

- Post-traumatic and preventative care
- Anatomically shaped stirrups ensure a perfect fit
- Unique strapping system prevents the forward shift of the talus, protecting the anterior talo-fibular ligament
- Low profile uprights easily fit under sport or dress shoes

Item sold as one (1) brace.

Indications

- Distorsion trauma of the upper ankle joint
- Insufficiency of the capsules and ligaments of the upper ankle joint
- Prevention of chronic instability

Suggested Code: L1902

Sizes + Product Codes

Size	Product Code	Women	Men
Left Brace			
S	BAS1L	5 - 8	< 7
M	BAS2L	9 - 12	7.5 - 10.5
L	BAS3L	13+	11 - 14

Size	Product Code	Women	Men
Right Brace			
S	BAS1R	5 - 8	< 7
M	BAS2R	9 - 12	7.5 - 10.5
L	BAS3R	13+	11 - 14



Elements of stability
U-shaped heel plate is anatomically shaped to the medial/lateral sides of the foot.