



Body Armor[®]

Cast Shoe

The redesigned Body Armor[®] Cast Shoe features improvements in function and design. It is specifically engineered to shield the cast from moisture, shock and everyday wear. The shoe can be used with all casting materials and bulky bandages while sporting a fresh, modern design that will appeal to all patients.

www.darcointernational.com





Bilateral



Light weight

Indications

- Cast protection
- Bulky bandages protection

Contraindications

- Open fractures

Suggested Code: L3260

The Benefits

Intentionally Designed Features

The design features sturdier low profile buttons, a more secure locking mechanism, a stronger bungee system, and a redesigned body that better adapts to the foot.

Built-in Protection

EVA Unibody Design provides unparalleled shock absorption while protecting the cast from the elements. Built in toe guard protects the toes and prevents debris from entering the toe area.

Ease of Use

Unique bungee closure system is easy to use, easy to adjust and effectively secures the cast in the shoe.



Sizes + Product Codes

Size	Insole				SKU #
	A	B	C	D	
XS	9.2	3.9	3.5	13.5	BACS0
S	10	4.4	4.0	14	BACS1
M	10.6	4.6	4.2	16	BACS2
L	11.4	5.1	4.6	18.5	BACS3
XL	12.6	5.5	4.8	20	BACS4
XXL	13.7	5.9	5.1	23	BACS5

