

CBO™

Open Toe Cast Boot

The CBO™ is a low profile multi use cast boot intended to protect casts of all types and heavy compression bandages.



Fits right or left foot



The Benefits

- **Heavy duty canvas** upper that protects the cast from the elements.
- Multilayer **EVA sole provides extra cushioning** for ease of ambulation and durability.
- The **rocker sole helps ease pressure** when ambulating and allows an easy roll through the gait pattern.
- Durable **hook and loop fastener straps seat and secure the casted foot** and keep the boot snug around the cast/bandaging.
- The CBO™ is **lightweight and low profile**.

Item sold as one (1) boot.

Suggested Code: L3260

Indications

- Fracture and cast shoe
- Bulky bandages

Sizes + Product Codes

Size	Children	Women	Men	Product Code
Pediatric	12 - 1	-	-	CBOP
XS	-	< 5	-	CBO0
S	-	5.5 - 7.5	-	CBO1
M	-	8 - 10	6 - 8.5	CBO2
L		-	9 - 11	CBO3
XL		-	11.5 - 14	CBO4