

# GTSTM

## Great Toe Splint

An adjustable, postoperative splint. Indicated for relief of symptoms associated with hallux valgus and other great toe joint conditions. Manufactured with lightweight neoprene and nylon, the GTSTM can be worn day and night, and its form-fitting design easily accommodates suitable footwear.



### Indications

- Hallux valgus and other great toe joint conditions

### The Benefits

- Maintains alignment
- 24 hour support
- Helps relieve painful bunions
- Can be worn postoperatively

*Item sold as one (1) splint.  
Available in right & left options.*

### Sizes + Product Codes

Size	Avg. Shoe Size	Product Code	
		Right	Left
Small	W 5 - 7.5 / M 3.5-6	GTS1R	GTS1L
Large	W 8 - 11 / M 6.5-9.5	GTS2R	GTS2L

**Suggested Code: L3100**



Fits comfortably under sock to accommodate most footwear

Comfortably corrects and maintains toe alignment