



The Body Armor® Heel Reliever is a versatile heel suspension boot designed to deliver targeted offloading and effective pressure relief throughout all phases of wound care—from prevention to the treatment of chronic pressure heel ulcers. Created for patients on bedrest or in a supine position for extended periods, it helps prevent decubitus heel ulcers (also known as pressure ulcers or bedsores) by redistributing pressure from the heel to the calf.

DARCO



New

Body Armor® Heel Reliever
Heel Protector



Bilateral



For diabetics



Time saving

Indications

- For the prevention and healing of pressure ulcers on the heel, in diabetic foot syndrome or decubitus
- Score of 14 or less on the Braden scale
- Diabetes with neurological sensory deficits
- Absence of palpable pulse on the foot
- Paraplegia
- Prolonged immobility or coma

Contraindications

- Wounds on the calf
- Spasticity
- Contraction of the ankle

Suggested Code: E0191

The Benefits

Targeted Relief + Offloading

- Offloads the heel by redistributing pressure to the calf
- Extra foam pads for customizable offloading




Uncompromised Protection

- Built-in positioning wedge prevents hip rotation
- Extended calf pad protects the knee from hyperextending



Thoughtfully Designed to Promote Healing

- Customizable design and accessories allow for a personalized fit
- Seamless inner surface reduces friction and shear
- Beveled edges to reduce contact pressure
- Lateral ventilation for improved air circulation
- Machine washable without losing shape or support



Color: Gray-Blue

Can be worn on the right or left foot

Size	Shoe Size	Calf Circumference (in)	Product Code
Universal	Men's 14 / Women's 15	9.8 - 21.3	BAHR



Version 2-BAHR_04_2026

DARCO International Inc.

810 Memorial Boulevard Huntington, WV 25701 | Toll Free: 800.999.8866



www.darcointernational.com

DARCO