



Night Splint Wedge

Posterior Night Splint

Designed to alleviate heel pain associated with plantar fasciitis and achilles tendonitis. It features two 1/2 inch EVA wedges, each with 3° of inclination, that can be removed or added to adjust the desired dorsiflexion of the foot. The natural cotton terry liner absorbs moisture and ensures a cool, comfortable fit and a non-skid sole for minimal ambulation.



Indications

- Plantar fasciitis
- Heel spurs
- Achilles tendon injuries
- Ankle contracture

The Benefits

- **Removable wedges** are included to create varying degrees of angle within the foot bed.
- **Fast and easy** to use with a hook and loop fastener.
- **Cotton terry liner** is soft, comfortable and absorbent.
- **Can be worn on the right and left.**
- Provides a **sustained stretch to the plantar fascia.**

Sizes + Product Codes

Size	Women	Men	Product Code
Small	Up to 6.5	-	NW1B
Medium	7 - 10	5 - 8	NW2B
Large	10.5 - 13	8.5 - 11	NW3B
X-Large	13.5+	11.5 - 14	NW4B



Item sold as one (1) splint. Suggested Code: L4396