



New

Body Armor® PFA

Partial Foot Orthosis

The Body Armor® PFA Walker is a lightweight, protective orthosis designed for post-surgical healing after partial foot amputation. Featuring a rigid outer shell with a soft, air-pump lining for even pressure and comfort, it's easy to apply with a wide opening and available in three sizes (M–XL).

The PowerPod EVA outsole makes it 25% lighter than previous models while delivering traction, shock absorption, and slip resistance. A ventilated honeycomb design and beveled construction encourage airflow, natural posture, and mobility during recovery.

DARCO

The Benefits

25% Lighter

- Innovative materials and strategic components make it 25% lighter than previous models

Accommodating Fit

- 3 size options (M-XL) for an optimized fit
- wide opening to accommodate swelling and for easy on-and-off
- Integrated Air Pump for even pressure distribution

PowerPod Outsole

- High-traction, EVA outsole pods grip securely for shock absorption and anti-slip
- Promotes natural posture and ergonomic movement
- Prevents delamination to improve product longevity



Wide Opening for easy on-and-off



Integrated Air Pump for even pressure distribution



Bilateral



For diabetics



Lightweight

Indications

- Postoperatively for partial foot amputations (e.g. Transmetatarsal, Chopart or Lisfranc amputations) due to:


- > Diabetic foot syndrome,
- > PAOD (peripheral arterial occlusive disease)
- > Trauma

Contraindications

- No Contraindications are known

Suggested Codes: L4360/ L4361

Adult Sizes + Product Codes



Size	Residual Limb Length (in)	Insole in inches			Product Code
		B	C	E	
M	7.7	5.2	3.7	7.9	BAPFA2
L	8.5	5.7	4.1	8.7	BAPFA3
XL	9.1	5.9	4.5	9.3	BAPFA4



PowerPod Outsole for stability and a secure fit