



ProMotion® Plus Orthotic Insoles


Designed for active individuals suffering from plantar fasciitis and other painful foot conditions, the ProMotion® Plus insoles provide long-lasting shock absorption, arch support, and heel padding.

Suggested Code: L3040



The Benefits

- **Full length PORON® cushioning foam** layer provides superior shock absorption and will not flatten or bottom out over time.
- **Smooth, anti-microbial polyester top layer** allows for freedom of movement, prevents sock bunch-up, and controls odor.
- **EVA covered polypropylene heel and arch support.** Provides a firm support that cradles the heel and arch area. The EVA surface prevents the insole from shifting and promotes its durability.



Indications

- Over-pronation
- Plantar fasciitis
- Heel pain
- Foot and arch pain
- Ankle pain
- Knee pain

Sizes + Product Codes

Item sold as pair of insoles.

Size	3	4	5	6	7	8	9	10	11	12	14
Men's (US)	-	-	5 - 5.5	6 - 6.5	7 - 7.5	8 - 8.5	9 - 9.5	10 - 10.5	11 - 11.5	12 - 13.5	14 - 14.5
Women's (US)	5 - 5.5	6 - 6.5	7 - 7.5	8 - 8.5	9 - 9.5	10 - 10.5	11 - 11.5	12 - 12.5	13 - 13.5	-	-
Product Code	PMP03	PMP04	PMP05	PMP06	PMP07	PMP08	PMP09	PMP10	PMP11	PMP12	PMP14