



Softie™

Post-Op Shoe

The DARCO Softie™ Post-Op Shoe is a cost-effective product that provides rigidity and control of the forefoot after soft tissue procedures or trauma to the foot.



The Benefits

- **Tri-laminated upper** is made of a lightweight, breathable poly-cotton blend that conforms to bony abnormalities
- **Semi-rigid metatarsal shank** provides postoperative healing and better protection after osteotomies.
- **Easy closure- two hook and loop fastener straps** allow for simple entry/exit and fastens anywhere along the surface of the soft upper for even compression and security.



Indications

- Postoperative for stabilization and off-loading
- Postoperative healing
- Trauma care
- Presence of K-Wires

Item sold as one (1) shoe. Fits right or left foot.

Suggested Code: L3260

Sizes + Product Codes

Size:	Product Code			
	Women	Men	Women	Men
Small	4 - 6	6 - 8	STW1B	STM1B
Medium	6.5 - 8	8.5 - 10	STW2B	STM2B
Large	8.5 - 10	10.5 - 12	STW3B	STM3B
X-Large	-	12.5 - 14	-	STM4B